



This document contains the confidential and proprietary information of Nanotronics Imaging, Inc. Neither this document nor any of the information contained is to be reproduced, distributed, used or disclosed, either in whole or in part, except as specifically authorized by Nanotronics Imaging, Inc.

nSpec Version 0.22.1.1

Release Date: 11 Nov 2022

Documentation Updated: 07 Nov 2022


Major Features: Interlock Bug Fixes, Custom Exporter (JSON, KLARF, XML)

- [Overview](#)
- [Upgrading to v0.22.1.1](#)
- [New Features](#)
 - [Highlights](#)
 - [NSPEC-5594: New Custom Exporter Outputs \(KLARF, JSON, XML\)](#)
 - [New Features Changelog](#)
- [Bug Fixes](#)
 - [Highlights](#)
 - [NSPEC-6267: Interlocks don't re-engage Stage](#)
 - [Changelog](#)
- [Appendices](#)
 - [Custom Exporter](#)
 - [Overview](#)
 - [Core Features](#)
 - [Embedding math and custom calculations](#)
 - [Supported Exports](#)
 - [GUI parameters](#)
 - [Referencing Prior Analyses](#)
 - [Configuring the nJson File](#)


Overview

This patch build contains two main highlights:

1. Fixes for some Interlock Systems
2. Additional export file options for the Custom Exporter

 As a patch build, QA testing has been limited to new feature tests, and did not include all regression tests.

Upgrading to v0.22.1.1

 Library Update Not Required

If upgrading from a version more than 1 release prior, please reference all intermediate release notes for upgrade steps *for each version*.

New Features

Highlights

NSPEC-5594: New Custom Exporter Outputs (KLARF, JSON, XML)

The Custom Exporter launched in nSpec version v0.22.1.0 with support for CSV exports. With this change the Custom Exporter now supports KLARF, JSON, and XML exports as well. For full documentation on this feature please see the page at [Custom Exporter](#)

New Features Changelog

T	Key	Release Notes Summary
	NSPEC-5594	Add KLARF support for the Custom Exporter
	NSPEC-6058	Autoloader Prealignment Speed Optimizations
	NSPEC-6180	Add JSON support for the Custom Exporter
	NSPEC-6181	Add XML support for the Custom Exporter

[4 issues](#)

Bug Fixes

Highlights

NSPEC-6267: Interlocks don't re-engage Stage

When the interlock on an enclosure unlocks, power is killed to the stage for safety reasons. However, when the door is shut and the interlock is locked again, the stage is not re-engaged until after the first motion. This causes a few problems:



- The first scan after the interlock was touched was captured at the incorrect position
- If the first motion is for an alignment job, first alignment mark will be missed

These problems have now been fixed! After open an interlock-equipped system door and closing it, stage motion will be immediately reset.

 Make sure to set the program option in the PMD category called `Use Interlock Recovery` to 1

Changelog

T	Key	Release Notes Summary	Affected Releases
---	-----	-----------------------	-------------------

-  [NSPEC-6267](#) Interlocks Don't Re-Engage the Stage 0.22.0.0, 0.22.1.0
-  [NSPEC-6292](#) Auto light adjust does not correct to set intensity value after first scan (cassette scanning) 0.22.0.0, 0.22.1.0

[2 issues](#)

Appendices

Custom Exporter

- [Overview](#)
- [Core Features](#)
 - [Embedding math and custom calculations](#)
 - [Supported Exports](#)
- [GUI parameters](#)
- [Referencing Prior Analyses](#)
- [Configuring the nJson File](#)

Overview


The Custom Exporter allows users to build highly customizable exports that report on all the desired defect information collected during prior analysis. Like the custom reporter and Dislocation Analyzer, the Custom Exporter also requires an nJson file loaded with relevant SQL scripts.

Core Features

CSV column headers are now completely configurable - no matter what nSpec calls a defect attribute, we have the ability to rename them to align with the customer vernacular. For example, within nSpec databases we might refer to defect width as “W”, but in this export we can rename that to “width” or “defectWidth”, or anything else that is required for smooth integration with the customer application.

Embedding math and custom calculations

In addition to simply exporting columns directly, the custom exporter will allow us to perform any operation that is legal within SQL code. This gives us enormous flexibility for reporting exactly what customers need.

 For more information on embedded calculations via the custom exporter, please contact the applications team

Supported Exports

Export File	Expected Release
CSV	v0.22.1.0
JSON	v0.22.1.1

XML	v0.22.1.1
KLARF (Device Inspection Scans)	v0.22.1.1
KLARF (All Scans)	TBD

GUI parameters




Analysis Parameters

Analysis: Custom Exporter

Group: Default

Description: Custom Exporter

Pos	Name	Type	Low	PV	High	Default	Picklist	Desc
01	Base Analysis IDs Mode	PickList		Relative		Relative	Absolute, Relative	the way next param (analysis ID) is treated
02	Base Analysis IDs	String						Coma Separated Analysis IDs
03	Export File Folder	PickFolder		\\4622105182\ProgramData\ConfigFiles\Analysis\Custom Exporter	...	\\462210518...		Choose the export file Folder
04	Custom Filename Suffix	String						Custom Filename Suffix
05	Export Setting	PickFile		\\4622105182\ProgramData\ConfigFiles\Analysis\Custom Exporter\Setting\Sample.njson	...	\\462210518...	njson	Custom Exporter File Path, if left empty it will use t

Param Name	Options	
Base Analysis IDs Mode	Absolute, Relative	<div>  See Post-Analysis Referencing section (Link) </div>
Base Analysis IDs		<div>  See Post-Analysis Referencing section (Link) </div>
Export File Folder		Select a folder to save the export file
Custom Filename Suffix	Optional	If non-empty, this suffix will be directly added to filename. If using nSpec v0.23.1.0 or earlier, the suffix will be added following an underscore.
Report Setting		Select a nJson file of the Report setting <div>  If you have a new request for a new functionality, please request support from support@nanotronics.co </div>

Referencing Prior Analyses

For guidance on how to set up post-analysis referencing, visit the appendix page for [Post-Analysis Referencing \[EXTERNAL\]](#)

Configuring the nJson File

The nJson file must be written and configured by Nanotronics Technical Service or Applications Engineers. To request a new custom reporting application that may require assistance, please contact support@nanotronics.co to request support.

For nJson formatting guidance, visit [nJson Format: Custom Exporter](#)